




CHAPTER 2 EXPERIMENTAL TECHNIQUES

11C UNDARGA



MEASUREMENT



VARIABLE	APPARATUS
• Time	• Stopwatch or Clock
• Temperature	• Thermomemeter (liquid in glass, thermistor or thermocouple)
• Mass	• Balance

**MEASURING VOLUME: BEAKER ,
BURETTE, PIPETTE , MEASURING
CYLINDER , GAS SYRINGE**




PURITY

PAPER CHROMATOGRAPHY IS USED TO SEPARATE MIXTURES OF SOLUBLE SUBSTANCES. THESE ARE OFTEN COLOURED SUBSTANCES SUCH AS FOOD COLOURINGS, INKS, DYES OR PLANT PIGMENTS. PAPER CHROMATOGRAPHY:

DROP SUBSTANCE TO CENTER OF FILTER PAPER AND ALLOW IT TO DRY
DROP WATER ON SUBSTANCE, ONE DROP AT A TIME



PAPER + RINGS = CHROMATOGRAM.
SUBSTANCES TRAVEL ACROSS PAPER AT DIFFERENT RATES WHICH IS WHY THEY SEPARATE INTO RINGS
METHOD WORKS BECAUSE DIFFERENT SUBSTANCES TRAVEL AT DIFFERENT LEVELS OF ATTRACTION TO IT



- STATIONARY PHASE IS MATERIAL ON WHICH SEPARATION TAKES PLACE
- MOBILE PHASE CONSISTS OF THE MIXTURE YOU WANT TO SEPARATE, DISSOLVED IN A SOLVENT.

RF VALUES


R_f values



**CHROMATOGRAMS OF PURE AND IMPURE SUBSTANCES
A PAPER CHROMATOGRAM CAN BE USED TO
DISTINGUISH BETWEEN PURE AND IMPURE
SUBSTANCES:**

- **A PURE SUBSTANCE PRODUCES ONE SPOT ON THE CHROMATOGRAM**
- **AN IMPURE SUBSTANCE PRODUCES TWO OR MORE SPOTS**

Water of Crystallisation




01 FILTRATION

- MIXTURE GOES IN A FUNNEL WITH FILTER PAPER, INTO A FLASK.
- RESIDUE IS INSOLUBLE AND STAYS AT TOP.
- FILTRATE GOES THROUGH

02 CRYSTALLISATION

- SOME WATER IN THE SOLUTION IS EVAPORATED SO SOLUTION BECOMES MORE CONCENTRATED.
- A DROP IS PLACED ON A SLIDE TO CHECK IF CRYSTALS ARE FORMING.
- SOLUTION IS LEFT TO COOL AND CRYSTALLISE.
- CRYSTALS ARE FILTERED TO REMOVE SOLVENT.



Filtration process


03 FRACTIONAL DISTILLATION

- REMOVES A LIQUID FROM A MIXTURE OF LIQUIDS, BECAUSE LIQUIDS HAVE DIFFERENT B.P.S
- MIXTURE IS HEATED TO EVAPORATE SUBSTANCE WITH LOWEST B.P.
- SOME OF THE OTHER LIQUID(S) WILL EVAPORATE TOO.
- A MIXTURE OF GASES CONDENSE ON THE BEADS IN THE FRACTIONAL COLUMN.
- SO THE BEADS ARE HEATED TO THE BOILING POINT OF THE LOWEST SUBSTANCE, SO THAT SUBSTANCE BEING REMOVED CANNOT CONDENSE ON THE BEADS.
- THE OTHER SUBSTANCES CONTINUE TO CONDENSE AND WILL DRIP BACK INTO THE FLASK.

04 SIMPLE DISTILLATION

- IMPURE LIQUID IS HEATED
- IT BOILS, AND STEAM RISES INTO THE CONDENSER
- IMPURITIES ARE LEFT BEHIND.

CONDENSER IS COLD SO STEAM CONDENSES TO THE PURE LIQUID AND IT DROPS INTO THE BEAKER



METHODS

Separation Techniques

Method of Separation	What to Separate?
Magnetic Attraction	Magnetic materials from non-magnetic materials
Filtration	A solid from a liquid
Crystallization	A solid from a solution
Simple Distillation	A solvent from a solution
Fractional Distillation	Liquids from each other
Chromatography	Different substances from a solution
Solvent Extraction	A component from a mixture with a suitable solvent Water dissolves some salts, sugar White spirit dissolves paint Propanone dissolves grease, nail polish Ethanol dissolves glues, printing inks, scented substances

SOLVENT	IT DISSOLVES...
Water	Some salts, sugar
White spirit	Gloss paint
Propanone	Grease, nail polish
Ethanol	Glues, printing inks, scented substances